



**PA-02** Package-Port<sup>®</sup>  
**PA-03** Package Entry System

**Enables a leak tight entry for seal strength and leak testing of flexible packages**


For use with the TM Electronics BT-1000 Automated Package Tester and Fixtures, the TME Package-Port<sup>®</sup> system allows consistent penetration of flexible packages for seal strength and leak integrity testing. The flexible Package Port is attached to the package using a double-sided adhesive disk to create a leak tight path between the package and the port. This allows the operator to simply and repeatably insert the test probe into the package for restrained and unrestrained package testing

***Applications in medical, food and pharmaceutical package testing***

***Low cost and dependable, easy to use, decreases package testing time***

***Increases the reliability of the results and brings real time test data to the production line***

***Complies with ISO-11607 standards and ASTM test methods***



**TMElectronics, Inc.**

45 Main Street, Boylston, MA 01505 USA  
800-370-0501 or 508-869-6400