



PA-02 Package Ports  
PA-03 Package Port Adhesive

## PA-02 Package-Port® PA-03 Package Entry System

**Enables a leak tight entry for  
seal strength and leak testing  
of flexible packages**

For use with the TME BT *Integra-Pack* Package Test Instrument and the TM Electronics BT-1000 Automated Package Tester and Fixtures, the TME Package-Port® system allows consistent penetration of flexible packages for seal strength and leak integrity testing. The flexible Package Port is attached to the package using a double-sided adhesive disk to create a leak tight path between the package and the port. This allows the operator to simply and repeatably insert the test probe into the package for restrained and unrestrained package testing

***Applications in medical, food and pharmaceutical package testing***

***Low cost and dependable, easy to use, decreases package testing time***

***Increases the reliability of the results and brings real time test data to the production line***

***Complies with ISO-11607 standards and ASTM test methods***



**TM Electronics, Inc.**

**45 Main Street, Boylston, MA 01505 USA  
800-370-0501 or 508-869-6400**

**sales@tmelectronics.com**

**www.tmelectronics.com**